

# VitaScan LT - C6 Model

To learn more about the C6 contact our Sales Department at (800) 757-8450, option 1 or sales@vitacon.us

Item# 100525C6 Includes (as shown): Bladder Scanner Probe, 10" Color Antimicrobial Medical Grade Tablet Rolling cart with basket and holders

### Real-time Bladder Image

Software identifies and tracks bladder location in real-time
Visualize physiological anomalies that may affect volume calculation
Option to manually adjust measurement based on real-time image
Post-test quality indicator
User-selectable scan depth

## **Fully Automated**

Multi-plane scan from one position on the body (3D) Automatically measures volume and Post Void Retention (PVR) Dynamic volume display (optional setting)

# **Seamless Data Capture**

Prints full-page reports to standard office printers Real-time save report as PDF, JPEG or DICOM Standard file formats for EMR storing Real-time save locally or to Network

# **Training and Support**

Free, unlimited training sessions for life
Calibration not required
Five Year Warranty with Accidental Damage Coverage<sup>1</sup>
Service and support team located in Minnesota, USA



#### **TECHNICAL SPECIFICATIONS**

SCANNER PROBE v2			
Volume Range	0 - 1000mL	Detection Depth 100, 160	or 230mm
Accuracy <sup>2</sup>	±7.5% when > 100mL	Probe Dimensions 7.3 x 1.8 >	1.8 inches
Probe Frequency	3.40 MHz	Probe Weight (incl. cable)	0.79 lb.
IP rating	IPX7	Power Consumption	0.5A@5V
Scanning Method	Sector, 180 degrees	Operating Conditions +50°F	to +113°F
Rotation Positions 6, 12 or 24 rotating positions		10% - 80% relative humidity	
Sweep Angle	130 degrees	Storage Temperature22°F	to +122°F
MEDICAL GRADE TABLET			
Operating System	Win 10 Pro	Printer Standard Offi	ce Printers
Display 10" Color Touchscreen		Housing Antimicrobial Medical Grade	

<sup>&</sup>lt;sup>1</sup> Accidental Damage Coverage requires deductible for replacement components.

<sup>&</sup>lt;sup>2</sup> Accuracy specification in accordance with stated instructions, scanning on Vitacon phantom using software version 4.5.4.14 and above.